## Replacing the Default Footer

Last Modified on 27/03/2019 1:23 pm ACDT

This article is part of the systems administration guide. You will require administration access to view the pages mentioned in this article.

## Get the Link for the Image

Once the image is uploaded you will need to get the link for the image.

- 1. From the Images/Pictures screen, select the Show button next to the relevant image.
  - a. The image will display on screen.
- The link for the image will display in the web address bar and will be similar to http://cr.test.clintelsystems.com/system\_images/footer
- Note down everything after your IP or web address, in this case we ignore http://cr.test.clintelsystems.com and note down: system\_images/footer
- 4. If the image is too large, you can resize it in the browser by adding the following to the address bar.
  - a. To resize to 400px wide and 40px high i.ehttp://cr.test.clintelsystems.com/system\_images/footer/size\_with/400x40
  - b. You can play with different sizes and when you find one that works note for the link
    e.g. system\_images/footer/size\_with/400x40

## Edit the Printing Template to Include the Footer

Once the image is uploaded and your have the image link you can add the footer to a custom printing template. (There is no footer included by default). This process assumes you have already created a custom (XSL-FO) printing template. If you haven't, see Customising a XSL-FO Printing Template

- 1. Go to **Administration**.
- 2. Select Correspondence from the menu.
- 3. Select Printing Templates.
- 4. Select XSL-FO.

Template.
a. The Template document screen will display.
6. Select <b>Edit</b> button at the top of the screen.
a. The XSL-FO template code will display.
7. Find the following Snippet of code:
8. Add the link to the image that you located previously (Get the link for the image) to the following section of code:
9. Next, you need to ensure the footer area is big enough for your image.
10. Find this snippet of code and change the size of themargin-bottom attribute - highlighted below:
11. Select <b>Update Custom Template</b> button.
Note: Generate a letter in Correspondence against a Test/Dummy patient. Access a patient record,

creating a letter and then Printing it.. You may have to test via trial and error to get the image to size

correctly.

5. Under Custom Templates section, select**Show** button on the relevant Default Letter